

# Bullfrog Batch Reactors - Single

# Bullfrog Batch Reactors



The high pressure reactors designed for you.

Whatever reactions you want to run, our Bullfrog Batch Reactors are flexible and robust enough to meet your needs.

With a large variety of options and accessories, you'll be able to select a configuration that matches your purpose - whether you want to evaluate a different reaction every time, or carry out standard tests every day. And if our standard products don't provide exactly what you want, our unique customisation capability means that we can design a solution specifically for you.

Our entry-level reactor comes in a range of sizes from 10ml up to 1L, is constructed from 316 stainless steel and has operating ranges of up to 200C and 200bar. The 4-port, quick-release lid design allows you to add your own choice of peripherals, such as a sampling line, pressure gauge and bursting disc.



# Single Reactors



For users with more advanced needs, standard options include Hastelloy or titanium construction, a 7-port lid and PEEK, Teflon or glass liners. Of course we can also provide a range of solutions for stirring, temperature and pressure management, with full computer control and data logging if required.

Our customised reactors will be the answer when your application has higher temperature and pressure demands, or where you need integration with existing systems.

If you are interested in high throughput solutions, take a look at the information on our Multiple Bullfrog Batch Reactors

All our batch reactors are simple to use, have a compact footprint and represent excellent value for money. And, of course, our reactors are always fully safety tested and can be certified to meet relevant legislation.



